

(19) World Intellectual Property
Organization
International Bureau



(43) International Publication Date
22 April 2004 (22.04.2004)


PCT

(10) International Publication Number
WO 2004/033647 A2

- (51) International Patent Classification⁷: C12N 126 East Lincoln Avenue, Rahway, NJ 07065-0907 (US). BENNETT, Paul, B., Jr. [US/US]; 126 East Lincoln Avenue, Rahway, NJ 07065-0907 (US). COHEN, Charles, J. [US/US]; 126 East Lincoln Avenue, Rahway, NJ 07065-0907 (US).
- (21) International Application Number: PCT/US2003/031822
- (22) International Filing Date: 9 October 2003 (09.10.2003)
- (25) Filing Language: English
- (26) Publication Language: English
- (30) Priority Data: 60/418,017 10 October 2002 (10.10.2002) US
- (74) Common Representative: MERCK & CO., INC.; 126 East Lincoln Avenue, Rahway, NJ 07065-0907 (US).
- (81) Designated States (*national*): CA, US.
- (84) Designated States (*regional*): European patent (AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HU, IE, IT, LU, MC, NL, PT, RO, SE, SI, SK, TR).
- (71) Applicant (*for all designated States except US*): MERCK & CO., INC. [US/US]; 126 East Lincoln Avenue, Rahway, NJ 07065-0907 (US).
- (72) Inventors; and
- (75) Inventors/Applicants (*for US only*): XIA, Meng-hang [US/US]; 126 East Lincoln Avenue, Rahway, NJ 07065-0907 (US). CONNOLLY, Thomas, M. [US/US];
- Published:
— without international search report and to be republished upon receipt of that report
- For two-letter codes and other abbreviations, refer to the "Guidance Notes on Codes and Abbreviations" appearing at the beginning of each regular issue of the PCT Gazette.

(54) Title: ASSAY METHODS FOR STATE-DEPENDENT CALCIUM CHANNEL AGONISTS/ANTAGONISTS

DOSE-RESPONSE RELATIONSHIP FOR K⁺-STIMULATED CALCIUM INFLUX IN WILD TYPE HEK 293 CELLS AND CELLS STABLY TRANSFECTED WITH THE L-TYPE Δ 1C CHANNEL (C1-6-37-3)



(57) **Abstract:** Methods of identifying activators and inhibitors of voltage-gated ion channels are provided in which the methods employ cells transfected with a voltage-gated ion channel of interest and a corollary channel to control the membrane potential of the cells by changing extracellular ion concentration. This allows for more convenient, more precise experimental manipulation of these transitions, and, coupled with efficient methods of detecting the result of ion flux through the channels, provides methods that are especially suitable for high throughput screening.